

# Applications of magnetite nanoparticles for water treatment and for DNA and cell separation

Hai N.H., Chau N., Luong N.H., Anh N.T.V., Nghia P.T.

Center for Materials Science, Faculty of Physics, Hanoi University of Science, Hanoi, 334 Nguyen Trai, Hanoi, Viet Nam; Center for Life Science Research, Faculty of Biology, Hanoi University of Science, Hanoi, 334 Nguyen Trai Road, Hanoi, Viet Nam

**Abstract:** Magnetic nanoparticles with a diameter of 15 nm prepared by using the coprecipitation method have been applied to enforce the sedimentation of the solid waste, to adsorb the arsenic ions in water, to increase the DNA concentration by using a magnetic field for the electrochemical DNA sensor and to separate the helper CD4<sup>+</sup> T cells to determine the number of the cells in blood. A combination of magnetic nanoparticles and alum makes the solid waste in water under a magnetic field aggregate a dozen times faster than under the gravity alone. A concentration of 0.1 mg/l of arsenic in water was reduced to a value lower than the permissible concentration of 0.01 mg/l after few minutes of stirring. The particles functionalized with 3-aminopropyl triethoxysilane were used to enrich the DNA of the Herpes virus, which extended the sensitivity of an electrochemical sensor down to a concentration lower than nM/l. The particles coated with fluorescent-labeled antiCD4 antibody were used to count the helper CD4<sup>+</sup> T cells. The fluorescence signals of the particle/cell system were two times stronger than those of the fluorescence antiCD4 cell system. This can be used for the treatment of an HIV-infected patient with a simple fluorescent microscope.

**Author Keywords:** Arsenic removal; CD4<sup>+</sup> T cell; Cell separation; DNA separation; Magnetic nanoparticles

Year: 2008

Source title: Journal of the Korean Physical Society

Volume: 53

Issue: 3

Page : 1601-1606

Cited by: 2

Link: [Scopus Link](#)

Correspondence Address: Hai, N. H.; Center for Materials Science, Faculty of Physics, Hanoi University of Science, Hanoi, 334 Nguyen Trai, Hanoi, Viet Nam; email: [nhhai@vnu.edu.vn](mailto:nhhai@vnu.edu.vn)

ISSN: 3744884

Language of Original Document: English

Abbreviated Source Title: Journal of the Korean Physical Society

Document Type: Conference Paper

Source: Scopus

Authors with affiliations:

1. Hai, N.H., Center for Materials Science, Faculty of Physics, Hanoi University of Science, Hanoi, 334 Nguyen Trai, Hanoi, Viet Nam
2. Chau, N., Center for Materials Science, Faculty of Physics, Hanoi University of Science, Hanoi, 334 Nguyen Trai, Hanoi, Viet Nam

Nam

3. Luong, N.H., Center for Materials Science, Faculty of Physics, Hanoi University of Science, Hanoi, 334 Nguyen Trai, Hanoi, Viet Nam
4. Anh, N.T.V., Center for Life Science Research, Faculty of Biology, Hanoi University of Science, Hanoi, 334 Nguyen Trai Road, Hanoi, Viet Nam
5. Nghia, P.T., Center for Life Science Research, Faculty of Biology, Hanoi University of Science, Hanoi, 334 Nguyen Trai Road, Hanoi, Viet Nam

#### References:

1. Pankhurst, Q.A., Connolly, J., Jones, S.K., Dobson, J., (2003) J. Phys. D: Appl. Phys, 36, pp. R167
2. D. L. Leslie-Pelecky, V. D. Labhasetwar and R. H. Kraus, Jr., Nanobiomagnetics in Advanced Magnetic Nanostructures, edited by D. J. Sellmyer and R. S. Skomski (Kluwer, New York, 2005)Khuat, N.T., Nguyen, V.A.T., Phan, T.-N., Thach, C.V., Hai, N.H., Chau, N., (2008) J. Korean Phys. Soc, 52, p. 1323
3. Thach, C.V., Hai, N.H., Chau, N., (2008) J. Korean Phys. Soc, 52, p. 1332
4. Hai, N.H., Pirn, N.D., Luong, N.H., Chau, N., Chinh, H.D., Hoang, L.H., Leslie-Pelecky, D.L., (2008) J. Korean Phys. Soc, 52, p. 1327
5. Bruce, I.J., Sen, T., (2005) Langmuir, 21, p. 7029
6. Yean, S., Cong, L., Yavuz, C.T., Mayo, J.T., Yu, W.W., Kan, A.T., Colvin, V.L., Tomson, M.B., (2005) J. Mater. Res, 20, p. 3255
7. Ryan, K.J., Ray, C.G., (2004) Sherris Medical Microbiology, , 4th ed, McGraw Hill, New York
8. Tuan, M.A., Binh, N.H., Tarn, P.D., Chien, N.D., (2005) Comm. Phys, 15, p. 218
9. [http://www.scioncorp.com/pages/scion\\_image\\_windows.htm](http://www.scioncorp.com/pages/scion_image_windows.htm)